

EPA LIST N

EPA approved on List N, for use against SARS-CoV-2, the cause of COVID-19.

Offered By:

CLEANCIDE GERMICIDAL PRODUCTS

A SMARTER, SAFER WAY TO CLEAN + DISINFECT

Citric acid-based CleanCide germicidal wipes deliver powerful efficacy against pathogens — without the harmful, irritating effects of other chemical disinfectants — in an easy-to-use wipe format that enables one-step, no-rinse cleaning and disinfection.



POWERFUL DISINFECTION – BROAD EFFICACY

Proven, broad spectrum, hospital-grade antimicrobial efficacy — virucidal, germicidal, fungicidal and tuberculocidal. Meets OSHA requirements for bloodborne pathogens and EPA List N criteria for use against SARS-CoV-2, the virus that causes COVID-19.



CITRIC ACID-BASED — NON-IRRITATING

A better alternative to quat-, alcohol- or bleach-based disinfectants. EPA Category IV (lowest level of toxicity) requiring no precautionary statements or PPE — will not irritate skin, eyes or lungs. Highly biodegradable ingredients, leaving no harmful or irritating residues.



IDEAL FOR HARD, NON-POROUS SURFACES

Formulated for fast, effective cleaning and disinfection of hard, nonporous surfaces such as stainless steel, chrome, glass, plastic and vinyl. Ideal for use on door knobs, bed rails, monitoring equipment and other digital monitors, telephones, light switches, faucet handles and other high-touch surfaces.



ONE-STEP, NO-RINSE CLEANING & DISINFECTION

Proven efficacy in the presence of 5% organic soil load. 5-minute contact time for disinfection. Enables one-step cleaning and disinfection — no pre-cleaning or post-rinsing steps required.



HOSPITAL GRADE



NO QUATS OR BLEACH



MULTI-SURFACE

PRODUCTS

PRODUCT #	DESCRIPTION	PACKAGING
2130-02	CleanCide Ready-To-Use*	12 Quarts
2130-07	CleanCide Ready-To-Use*	5 Gal Pail
3130LG-12	CleanCide Wipes – 60 Count Canister	12 Canisters
3130C-160-CS	CleanCide Wipes – 160 Count Canister	12 Canisters
3130FP-80	80 count flow pack	each
3130B-400	CleanCide Wipes – 400 Count Tub	1 Tub



* EPA Designed for the Environment (DfE) product

† Sanitizing claim only applies to RTU liquid

CLEANCIDE EFFICACY DATA

EPA REGISTRATION NUMBER 34810-35

DISINFECTING

BACTERICIDAL

Bactericidal activity (in the presence of 5% blood serum)

3 MINUTES	ATCC
Carbapenem-resistant Klebsiella pneumoniae	ATCC BAA-1705
Enterococcus faecalis Vancomycin Resistant VRE	ATCC 51575
Escherichia coli O157:H7	ATCC 43888
Listeria monocytogenes	ATCC 19117
Pseudomonas aeruginosa	ATCC 15442
Salmonella enterica (formerly Salmonella Choleraesuis)	ATCC 10708
Staphylococcus aureus	ATCC 6538
Staphylococcus aureus - MRSA	ATCC 33592
Staphylococcus epidermidis - MRSE	ATCC 51625

VIRUCIDAL

In the presence of 5 % blood serum, is effective against

3 MINUTES	STRAIN AND/OR ATCC
Hepatitis B (HBV), Duck Hepatitis B (DHBV) as surrogate	Grimaud Strain
Hepatitis C (HCV), Bovine Viral Diarrhea Virus (BVDV) as surrogate	Strain NADL
Rhinovirus	Strain 1059
5 MINUTES	
Adenovirus Type 2	ATCC VR-846
Feline Calicivirus (surrogate for Norovirus)	Strain F-9
Herpes Simplex Virus Type 1	ATCC VR-733
Herpes Simplex Virus Type 2	ATCC VR-734
Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)	Strain HTLV-III B
Influenza A Virus	Strain Hong Kong H3N2
Influenza A Virus (formerly called swine flu)	Strain Pandemic 2009 H1N1
Respiratory Syncytial Virus	ATCC VR-26
Rotavirus	Strain WA, Univ. of Ottawa
Vaccinia Virus	ATCC VR-119

FUNGICIDAL

Fungicidal activity (in the presence of 5% blood serum)

5 MINUTES	ATCC
Candida albicans	ATCC 10231
Trichophyton interdigitale (T. interdigitale)	ATCC 9533

TUBERCULOCIDAL

Tuberculocidal activity at 20°C/68°F on previously cleaned hard, non-porous, inanimate surfaces

5 MINUTES	STRAIN
Mycobacterium bovis (Mycobacterium tuberculosis) (M. bovis) (TB)	Strain BCG (Organon Teknika Corp.)

SANITIZING*

60 SECONDS	ATCC
Klebsiella pneumoniae	ATCC 4352
FOOD CONTACT SANITIZER: 60 SECONDS	ATCC
Staphylococcus aureus	ATCC 6538
Escherichia coli	ATCC 11229

* Efficacy claims apply to RTU product only

CONTACT YOUR CLEANCLIDE REPRESENTATIVE TODAY!



FACILITY PRODUCTS
COMBINING NATURE + SCIENCE + INNOVATION



WEXFORDLABS

© 2020 Wexford Labs Inc. All Rights Reserved.

3021 W. FM 1093, Wallis, TX 77485
979.478.2489 ofc.

PO Box 1030, Wallis, TX 77485

sales@facilityproducts-sf.com

713.960.8554 fax